












Archive Explorer Toolbar

Icon	Shortcut Keys	Description
		Load the selected study or worklist.
	Ctrl + E	Explore the selected item by navigating down one level.
	Ctrl + Backspace	Navigate up one level.
	Home	Go to your home folder.
		Copy the selected data.
		Burn the selected data to CD.
		Delete the selected data.
		Stop loading the study in progress.
	F5	Refresh the Archive Explorer.
Sub-Select...		Specify the criteria for selecting items on the current list: every 2, every 3, every 5, every 10, or select all items, or none (deselect all items).
		Show patient reports.
		Show Sticky Notes.