

Viewer Shortcut Keys

Shortcut Keys	Command	Shortcut Keys	Command
P	Pan	Ctrl + Q	Open and run Cine/Stop and close Cine
Z	Zoom	Ctrl + R	Reports
Ctrl + 1	Confirm and load next study	Ctrl + S	Stack
Ctrl + 2	Confirm and open Archive Explorer	Ctrl + T	Teaching file
Ctrl + 3	Dismiss and load next study	Ctrl + U	Unlink all
Ctrl + 4	Dismiss and open Archive Explorer	Ctrl + V	Paste (images/graphics)
Ctrl + A	Automatic windowing	Ctrl + W	Swivel
Ctrl + B	Shutter	Ctrl + Y	Delete all graphics
Ctrl + C	Copy to clipboard (images/graphics)	Ctrl + Space	Key image preview
Ctrl + D	Dictate	Ctrl + Tab	Next viewbox
Ctrl + Shift + D	Multiple Dictation	Ctrl + Shift + Tab	Previous viewbox
Ctrl + E	Enlarge group	Alt + Space	Mark key image
Ctrl + F	Reset all	Alt + D	Reapply original display protocol
Ctrl + G	Hide all graphics	Alt + E	Enhance --> Sharp
Ctrl + H	Hide annotations	Alt + G	Hide system graphics (in MPR)
Ctrl + I	Pointer	Alt + I	Inverse
Ctrl + K	Sticky Notes	Alt + L	Line measurement tool
Ctrl + L	Automatic Link	Alt + P	Save presentation
Ctrl + M	Magnify	Alt + S	New session
Ctrl (while magnifying)	Adjust magnification factor inside magnifying glass	Alt + X	Dictate
Alt (while magnifying)	Adjust size of magnifying glass	Alt + F2	Exit
Shift (while magnifying)	Adjust window/level inside magnifying glass	Alt + 0	Open Patient Mini-Archive next to cursor position
Ctrl + N	Explore Archive	F1 - F8	Window/Level preset/on
Ctrl + O	Open Patient Mini-Archive	F11	Full screen
		F12	Film single image