

SENSITITRE™ GRAM POSITIVE PLATE FORMAT

Plate Code: **GPALL3F**

	1	2	3	4	5	6	7	8	9	10	11	12
A	CPT 0.12	CPT 0.25	CPT 0.5	CPT 1	CPT 2	CPT 4	AMP 0.25	AMP 0.5	AMP 1	AMP 2	AMP 4	AMP 8
B	CHL 2	CHL 4	CHL 8	CHL 16	AXO 8	AXO 16	AXO 32	AXO 64	DAP 0.5	DAP 1	DAP 2	DAP 4
C	CLI 0.5	CLI 1	CLI 2	GEN 4	GEN 8	LEVO 0.25	LEVO 0.5	LEVO 1	LEVO 2	LEVO 4	DT1	DT2
D	LZD 1	LZD 2	LZD 4	LZD 8	VAN 0.25	VAN 0.5	VAN 1	VAN 2	VAN 4	VAN 8	VAN 16	FOXS 6
E	PEN 0.12	PEN 0.25	PEN 0.5	PEN 1	PEN 2	PEN 4	PEN 8	TLA 0.015	TLA 0.03	TLA 0.06	TLA 0.12	TLA 0.25
F	MXF 0.25	MXF 0.5	MXF 1	MXF 2	MXF 4	OXA+ 0.25	OXA+ 0.5	OXA+ 1	OXA+ 2	SXT 2/38	SXT 4/76	POS
G	RIF 1	RIF 2	TGC 0.03	TGC 0.06	TGC 0.12	TGC 0.25	TGC 0.5	NIT 32	NIT 64	GEN 500	STR 1000	POS
H	TET 2	TET 4	TET 8	ERY 0.5	ERY 1	ERY 2	ERY 4	CIP 1	CIP 2	CIP 4	NEG	POS

ANTIMICROBICS

CPT	Ceftaroline
CHL	Chloramphenicol
CLI	Clindamycin
LZD	Linezolid
PEN	Penicillin
MXF	Moxifloxacin
RIF	Rifampin
TET	Tetracycline
TGC	Tigecycline
GEN	Gentamicin
ERY	Erythromycin
AXO	Ceftriaxone
VAN	Vancomycin
LEVO	Levofloxacin
OXA+	Oxacillin+2%NaCl
AMP	Ampicillin
TLA	Telavancin w/Tween mimic
NIT	Nitrofurantoin
CIP	Ciprofloxacin
DAP	Daptomycin
SXT	Trimethoprim / sulfamethoxazole
DT1	D Test 1
STR	Streptomycin
NEG	Negative Control
DT2	D Test 2
FOXS	Cefoxitin screen
POS	Positive Control