



How to Reconstitute Antibiotic Beyond Use Date: Amoxicillin

Step-By-Step: Amoxicillin 50mg/ml 80ml

1. Shake bottle to loosen powder
2. Use syringe to measure 59ml of water and add into Amoxicillin bottle.
Powder may clump – be sure to reconstitute all powder well
3. Insert bottle adaptor (if applicable) and place cap on bottle
4. Shake vigorously
5. Use Amoxicillin within 14 days. Store at either room temperature or preferred refrigeration (but not required)
6. Always shake well prior to administration



Scan the QR code or visit <https://go.illinois.edu/AntibioticReconstitutionVetMed> for a video demonstration of this process

Large Animal Clinic

Emergency, specialty, and primary care for horses and farm animals
1102 W. Hazelwood Drive, Urbana
(217) 333-2000

Small Animal Clinic

Emergency and specialty services for small companion animals
1008 W. Hazelwood Drive, Urbana
(217) 333-5300

South Clinic

Health and wellness care for dogs, cats, and zoological pets
2100 Goodwin Avenue, Urbana
(217) 244-2555

Midwest Equine

Equine-only practice that provides equine veterinary and farrier services
20844 Cumberland Road, Farmer City
(309) 928-2127